

In-House Lab Work Routing Slip



DUE DATE: _____

Seven (7) days is standard unless ok'd by Lab Assistant and/or Dr. Rodriguez

Patient Name: _____ Today's Date: _____

CIRCLE THE TYPE OF APPLIANCE BEING MADE:

Occlusal Splints

- Hard/Soft Splint (standard)
- Fully Rigid Splint
- Soft Splint

Whitening Trays

(upper lower both)

Essix (Invisalign-type) Retainer

(upper lower both)

Athletic Mouthguard (upper only)

- Physical contact: (light medium heavy)
- Color: (clear white red blue black)

I have provided either:

- **unwrapped PVS impression(s), or**
- **alginate impression(s), wrapped in moist (not soaked!) paper towel, and sealed within plastic bag.**

Signed: _____

- Only bubble/void free impressions will be accepted for lab work. If you are unsure, retake!
- PVS is great at hard tissue detail, while alginate is best at capturing both hard & soft tissue details.
- PVS impressions must be allowed to de-gas for 12-24h before pouring up.
- Alginate impressions can be poured up immediately, or within 24h (if stored properly).
- For athletic mouthguards, use alginate only and be sure and capture entire maxillary vestibule space.

Special Notes:

In-House Lab Work Routing Slip



DUE DATE: _____

Seven (7) days is standard unless ok'd by Lab Assistant and/or Dr. Rodriguez

Patient Name: _____ Today's Date: _____

CIRCLE THE TYPE OF APPLIANCE BEING MADE:

Occlusal Splints

- Hard/Soft Splint (standard)
- Fully Rigid Splint
- Soft Splint

Whitening Trays

(upper lower both)

Essix (Invisalign-type) Retainer

(upper lower both)

Athletic Mouthguard (upper only)

- Physical contact: (light medium heavy)
- Color: (clear white red blue black)

I have provided either:

- **unwrapped PVS impression(s), or**
- **alginate impression(s), wrapped in moist (not soaked!) paper towel, and sealed within plastic bag.**

Signed: _____

- Only bubble/void free impressions will be accepted for lab work. If you are unsure, retake!
- PVS is great at hard tissue detail, while alginate is best at capturing both hard & soft tissue details.
- PVS impressions must be allowed to de-gas for 12-24h before pouring up.
- Alginate impressions can be poured up immediately, or within 24h (if stored properly).
- For athletic mouthguards, use alginate only and be sure and capture entire maxillary vestibule space.

Special Notes:
